

# OraGRAFT®

## Cortical Mineralized Powder



The OraGraft product line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects.

Brand	<ul style="list-style-type: none"> <li>OraGraft</li> </ul>
Clinical Applications	<ul style="list-style-type: none"> <li>Periodontal Defects</li> <li>Ridge Maintenance</li> <li>Extraction Site Preservation</li> <li>Sinus Grafting</li> <li>Craniofacial Reconstruction</li> </ul>
Benefits	<ul style="list-style-type: none"> <li><b>Highest Safety Level:</b> Sterile with an SAL of <math>10^{-6}</math></li> <li><b>Convenient Storage:</b> Ambient temperature with a 3 year shelf life</li> <li><b>Graft Incorporation:</b> Osteoconductive graft with dense architecture allows for healing through "creeping substitution"</li> <li><b>Versatile:</b> Available in 5 size configurations</li> </ul>
Alternative Products <i>(with surgeon approval)</i>	<ul style="list-style-type: none"> <li>Cortical/Cancellous Mix</li> <li>Cancellous Mineralized Particulate</li> </ul>

# OraGRAFT®

## Cortical Mineralized Powder



### CORTICAL MINERALIZED POWDER

250 - 1000 microns

Size	Freeze-Dried
0.25 cc	GCI/20
0.50 cc	GCI/10
0.70 cc	GCI/8
1.20 cc	GCI/4
2.50 cc	GC

